









EASYCAP


EEG Recording Caps and Related Products


EASYCAP GmbH
Steingrabenstrasse 14
DE-82211 Herrsching
Germany

Tel +49 (0) 8152 3722-24
Fax +49 (0) 8152 3722-29
www.easycap.de
info@easycap.de

actiCAP with 35 + 2 cap positions for Jos van Berkum, MPI Nijmegen, 1/08

-  Green holders: Label 1-32, hard-wired Ch1 – Ch32
-  Pink holders: Label 1-32, hard-wired Ch65 – Ch96
-  Blue holder: Label & hard-wired Ref
-  White holders: Label 1-32, hard-wired Ch33 – Ch64
-  Transp. holders: Label 1-32, hard-wired Ch97–Ch128
-  Black holder: Label & hard-wired Gnd



 Electrode with holder for placement outside cap on skin.









EASYCAP



EEG Recording Caps and Related Products

EASYCAP GmbH
Steingrabenstrasse 14
DE-82211 Herrsching
Germany

Tel +49 (0) 8152 3722-24
Fax +49 (0) 8152 3722-29
www.easycap.de
info@easycap.de

actiCAP with 60 + 2 cap positions for Jos van Berkum, MPI Nijmegen, 1/08

-  Green holders: Label 1-32, hard-wired Ch1 – Ch32
-  White holders: Label 1-32, hard-wired Ch33 – Ch64
-  Pink holders: Label 1-32, hard-wired Ch65 – Ch96
-  Transp. holders: Label 1-32, hard-wired Ch97–Ch128
-  Blue holder: Label & hard-wired Ref
-  Black holder: Label & hard-wired Gnd



Electrodes with holders for placement outside cap on skin.









EASYCAP


EEG Recording Caps and Related Products

EASYCAP GmbH
Steingrabenstrasse 14
DE-82211 Herrsching
Germany

Tel +49 (0) 8152 3722-24
Fax +49 (0) 8152 3722-29
www.easycap.de
info@easycap.de

actiCAP with 123 + 2 cap positions for Jos van Berkum, MPI Nijmegen, 1/08

-  Green holders: Label 1-32, hard-wired Ch1 – Ch32
-  White holders: Label 1-32, hard-wired Ch33 – Ch64
-  Pink holders: Label 1-32, hard-wired Ch65 – Ch96
-  Transp. holders: Label 1-32, hard-wired Ch97–Ch128
-  Blue holder: Label & hard-wired Ref
-  Black holder: Label & hard-wired Gnd



Electrodes with holders for placement outside cap on skin.